

## THE STERIKING COVER BAGS

are manufactured of tough, transparent plastic laminate. The plastic covers are developed to protect sterile packs from dust, moisture and air-borne contamination. You can cover your packs in any clean, appropriate work area.



## THE STERIKING COVER BAGS

are to be used after sterilization. A sterilized pack is put into the Steriking Cover Bag, the heat-seal feature of which assures fast, reliable and secure closure. The recommended sealing temperature is 130—160°C.

The Cover Bags are peelable, facilitating aseptic opening of a pack. The simple and effective opening feature allows the contents to be quickly and easily removed.

## THE STERIKING COVER BAGS

extend the shelf-life of the sterilized items by providing an effective bacteria barrier.

The following standard sizes are available ex stock:

150 × 270 mm,	1 000 pcs/box	200 × 300 mm,	1 000 pcs/box
250 × 400 mm,	1 000 pcs/box	300 × 500 mm,	500 pcs/box
400 × 550 mm,	500 pcs/box	420 × 620 mm,	500 pcs/box



WIPAK MEDICAL  
P.O. Box 45  
SF-15561 NASTOLA  
FINLAND  
Tel. 358/18 62011  
TLX 16181 wipak sf  
Tfax 358/18-620 300

WIPAK MEDICAL  
18, bis rue d'Anjou  
F-75008 PARIS  
FRANCE  
Tel. direct 33/1-42 65 96 11  
33/1-42 65 86 94  
Tlx 650 297 copap

WIPAK MEDICAL  
P.O. Box 1259  
Siemensring 116  
D-4156 WILLICH 1  
FEDERAL REPUBLIC  
OF GERMANY  
Tel. 49/02154-1085-87  
Tlx 8531900  
Tfax 49/02154-42135

WIPAK MEDICAL  
Schoenerweg 56  
NL 6222 NX MAASTRICHT  
THE NETHERLANDS  
Tel. 31/43-631 300  
Tlx 56946 wipak nl  
Tfax 31/43-631 563

GREENCROWN  
PACKAGING LTD.  
24, Heatfield  
Stacey Bushes  
MILTON KEYNES MK 12 6HP  
GREAT BRITAIN  
Tel. 44/908-317 754  
Tlx 826714 grnwn g

WIPAK MEDICAL  
Upplagsvägen 27  
S-100-27 STOCKHOLM  
SWEDEN  
Tel. 46-8-7440375  
Tlx 11205 sand s