

Sterilization **Tray-Liners** for absorbing excessive moisture and for mechanical protections

Steriking® tray-liners are specially developed for a double function. They absorb condensation in surgical instrument trays during steam sterilization processes. Sterility of a processed pack can come to risk due to excessive wetness. The tray-liners eliminate effectively the lint problems that conventional linen sheets cause. Furthermore they protect instruments from direct contact with tray wires.

TL – Standard Tray-Liner

The tray-liner sheet is placed into a tray under the instrument. It replaces effectively the linen sheets and is a cost effective choice for regular use.

Order code	Size	Packing unit
TL25	25x30 cm	500 sheets
TL40	40x50 cm	500 sheets
TL50	30x50 cm	500 sheets



New!

TLWE – Tray-Liner With Edge

The tray-liner with edge protects the sterilized package against holes and damage, and it absorbs excessive condensation during autoclaving. It is placed outside of the instrument tray on the top of the innermost wrap. This strong protection sheet with high absorption capacity protects sharp instruments or tray edges from sticking through the wrapping material.

Order code	Inner size	Outer size with edges	Packing unit
TLWE5030	50x30 cm	60x40 cm	500 sheets
TLWE4825	48x25 cm	58x35 cm	500 sheets
TLWE2525	25x25 cm	35x35 cm	500 sheets



New!

TLN – Nonwoven Tray-Liner

The TLN tray-liners are highly recommended for use in combination with nonwoven and SMS type of wrapping materials. It is very soft and very strong and absorbs water 6 times its own weight. Due to the extreme absorption capacity the nonwoven tray-liners are an optimal solution with heavy instrument trays when condensation is very obvious.

Order code	Size	Packing unit
TLNW 1525	15x25 cm	3200 sheets
TLNW 3025	30x25 cm	1200 sheets
TLNW 3040	30x40 cm	1200 sheets
TLNW 3050	30x50 cm	600 sheets
TLNW 5040	50x40 cm	800 sheets



WIPAK MEDICAL

Finland: Wihuri Oy Wipak, P.O. Box 45, Wipaktie 2, FI-15561 Nastola
Tel. +358 3 468 311, Fax +358 3 468 3333
E-mail: info.finland@wipak.com

Netherlands: Wipak B.V., Nieuwsadterweg 17, NL-6136 KN Sittard
Tel. +31 46 420 2999, Fax +31 46 458 1311
E-mail: info.benelux@wipak.com

www.steriking.fi www.wipak.com